

# Environment of Care product portfolio



**Sani-Cloth<sup>®</sup> Prime**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **1** MINUTE

DISINFECTANT  
QUAT/  
ALCOHOL  
FORMULATION

KILLS **56** MICROORGANISMS



A 1 MINUTE OVERALL CONTACT TIME, INTERMEDIATE LEVEL DISINFECTANT FOR ALL YOUR HIGH-TOUCH, HOUSE-WIDE SURFACE DISINFECTION NEEDS

**Sani-HyPerCide<sup>®</sup>**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **1** MINUTE

EFFECTIVE AGAINST **C. DIFF**  
**5** MINUTES

KILLS **24** MICROORGANISMS



STANDARDIZE YOUR DISINFECTION PROTOCOLS WITH THIS POWERFUL H<sub>2</sub>O<sub>2</sub> FORMULA THAT KILLS *CLOSTRIDIODES DIFFICILE* WITHOUT COMPROMISING COMPATIBILITY

**Super Sani-Cloth<sup>®</sup>**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **2** MINUTES

#1 DISINFECTANT WIPE IN HEALTHCARE<sup>1</sup>

DISINFECTANT  
QUAT/  
ALCOHOL  
FORMULATION

KILLS **32** MICROORGANISMS



IDEAL FOR DAILY USE IN FAST-PACED ENVIRONMENTS THAT REQUIRE SHORT CONTACT TIMES AND BROAD COVERAGE OF MICROORGANISMS

**Sani-Cloth<sup>®</sup> AF3**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **3** MINUTES

DISINFECTANT  
QUAT  
ALCOHOL FREE  
FORMULATION

KILLS **45** MICROORGANISMS



PROTECTING YOUR PATIENTS HAS NEVER BEEN EASIER. FRAGRANCE FREE – IDEAL FOR USE AROUND PATIENTS AND STAFF, ESPECIALLY THOSE WITH RESPIRATORY SENSITIVITIES

**Sani-Cloth<sup>®</sup> Plus**  
GERMICIDAL DISPOSABLE CLOTH

DISINFECTS IN **3** MINUTES

DISINFECTANT<sup>2</sup>  
QUAT  
FORMULATION

KILLS **18** MICROORGANISMS



IDEAL FOR USE IN ALTERNATE CARE SETTINGS INCLUDING PHYSICIAN AND DENTAL OFFICES

**Sani-Cloth<sup>®</sup> Bleach**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **4** MINUTES

DISINFECTANT  
1:10 BLEACH  
DILUTION  
FORMULATION

KILLS **52** MICROORGANISMS



IDEAL FOR DISINFECTING AREAS CONTAMINATED WITH *CLOSTRIDIODES DIFFICILE* SPORES, *CANDIDA AURIS*, AND NOROVIRUS

**Easy Screen<sup>®</sup>**  
CLEANING WIPE

IDEAL FOR ANY EQUIPMENT  
REQUIRING 70% IPA

CLEANING  
70% ISOPROPYL  
ALCOHOL  
FORMULATION



A COMPATIBLE 70% IPA TOUCHSCREEN CLEANER THAT WORKS AT THE SPEED OF TECHNOLOGY<sup>3</sup>

<sup>1</sup>Kills SARS-CoV-2 (Coronavirus), the virus that causes COVID-19, on hard, nonporous surfaces excluding Easy Screen<sup>®</sup> wipes.

<sup>2</sup>Effective against *Candida auris* (*C. auris*).

<sup>3</sup>Market leadership claims based on Q1 2021 Clarivate Data; <sup>2</sup>Sani-Cloth<sup>®</sup> Plus Germicidal Disposable Wipe compares to competitor Quat/Low-Alcohol disinfectant products; <sup>3</sup>Cleaning the screen of your mobile device with a wipe containing alcohol may reduce the lifespan of any anti-fingerprint coating. Please check your device manufacturer's cleaning instructions before using Easy Screen<sup>®</sup> wipes and all other PDI disinfectant products. Easy Screen<sup>®</sup> wipes are to be used on energized electronic equipment only.

## Sani-Cloth® Prime

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P25372	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P24284	7.5" x 15"	70/can 6 cans/case
Large Individual Packets	H06182	5" x 8"	50/box 10 boxes/case
X-Large Individual Packets	U13195	11.5" x 11.75"	50/box 3 boxes/case

## Sani-HyPerCide®

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P27372	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P26584	7.5" x 15"	65/can 6 cans/case

## Super Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	Q55172	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P86984	7.5" x 15"	75/can 6 cans/case
Large Individual Packets	H04082	5" x 8"	50/box 10 boxes/case
X-Large Individual Packets	U87295	11.5" x 11.75"	50/box 3 boxes/case
Softpack	A22480	8.2" x 9.8"	80/pack 9 packs/case

## Sani-Cloth® AF3

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P13872	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P72584	7.5" x 15"	75/can 6 cans/case
Large Individual Packets	H59200	5" x 8"	50/box 10 boxes/case
X-Large Individual Packets	U27500	11.5" x 11.75"	50/box 3 boxes/case
Softpack	M8063S80	8.2" x 9.8"	80/pack 9 packs/case
Pail	P1450P	7.5" x 15"	160/pail 2 pails/case
Refill (for Pail)	P2450P	7.5" x 15"	160/refill 2 refills/case

## Sani-Cloth® Plus

GERMICIDAL DISPOSABLE CLOTH

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	Q89072	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	Q85084	7.5" x 15"	65/can 6 cans/case

## Sani-Cloth® Bleach

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Clinical Size Wipes Canister	P84172	6" x 5"	160/can 12 cans/case
Large Canister	P54072	6" x 10.5"	75/can 12 cans/case
Extra-Large Canister	P25784	7.5" x 15"	65/can 6 cans/case
Large Individual Packets	H58195	5" x 7"	40/box 10 boxes/case
X-Large Individual Packets	U26595	11.5" x 11.75"	40/box 3 boxes/case
Pail	P7007P	7.5" x 15"	160/pail 2 pails/case
Refill (for Pail)	P700RF	7.5" x 15"	160/refill 2 refills/case

## Easy Screen®

CLEANING WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Medium Canister	P03672	6" x 9"	70/can 12 cans/case
Packet	H01050	5" x 8"	50/box 10 boxes/case

### Educational and point-of-care accessories.

Developed to assist you in facility protocol compliance.



PDI Sani-Bracket®  
for wall/pole mounting of a canister



Sani-Canister Caddy®  
for placement on flat surfaces (bilingual)



Sani-TAG® Equipment ID System  
for identification of equipment to ensure the proper  
wipe is used on a particular piece of equipment



Educational Wall Charts and  
signage (bilingual)



Scan here to contact your local PDI representative.